

## How to Personalise your Tester Blotter Strip

- 1.** Select the quantity of Tester Blotter Strips you require from the drop down menu.
- 2.** Click **PERSONALISE**. You can now view the shape of the Tester Blotter Strip that you will be dropping your logo onto. Your logo can be in JPG, PNG or GIF format.
- 3.** Click anywhere on the greyed out Tester Blotter Strip Image and a '**Drop a file here**' box will appear. Either drag and drop your logo into this box or click on the 'choose a local file' button to select your logo. Once your logo has been selected scroll down within this box and you will see a blue button **ADD**. Press this button and your logo will attach to the Tester Blotter Strip.
- 4.** You will now see your logo within the Tester Blotter Strip. Position your logo by clicking on your logo and dragging to the position you prefer. Use the controls to rotate your logo / increase or decrease its size as you wish.
- 5.** Click **APPROVE**. Your personalised Tester Blotter Strip order will be assigned to your **SHOPPING BASKET**. Press **CANCEL** to re-start the process.